MARSXpress 2.0

Scientific Microwave System



The MARSXpress 2.0™ system brings simplicity, flexibility and accessibility to every laboratory. As the most cost-effective solution available, MARSXpress 2.0 with XpressIntelligence will take care of your sample preparation needs.

Simplicity

Effortless programming with XpressIntelligence. Use XpressStart methods to start a digestion: select the XpressStart temperature and press start; it's that easy. The system automatically determines the rest. The world's most popular vessel, MARSXpress, require no tools for the simple three-piece assembly and range in size from 20 mL to 130 mL.

Flexibility

Workflows that work for you. Batch from 1, 2, 3, or up to 40 samples. Whether digesting a single sample for rush needs, performing a quick screen for optimal conditions, or digesting 40 samples in a batch, simply load the vessels and press start.

Accessibility

Information like never before. Check on your batch status from across the lab, watch digestions take place through the largest window available, review run data, generate reports, and more. Information is at your fingertips.

Learn more at cem.com/marsxpress-2

CEM Advancing Science

cem.com

© 2026 CEM Corporation. All Rights Reserved.

EDGE

Automated Solvent Extraction System





The EDGE is faster than Soxhlet, more automated than QuEChERS, and simpler than other solvent extraction systems. Extract a wide range of samples and sizes at least 3 times faster than other pressurized fluid extractors, including filtering, cooling, and washing. It has revolutionized the extraction process for sample preparation.

Simplicity

EDGE combines the speed of pressurized fluid extraction with the simplicity of Q-Cup® technology. It is designed to make sample preparation straightforward: load your samples, select a method, and press Start. EDGE automates solvent addition, extraction, filtration, cooling and washing — delivering a particle-free extract ready for analysis.

Flexibility

EDGE supports a broad range of sample types and extraction needs, including:

- Environmental samples
- Food and agricultural products
- Plastics and consumer products
- Pesticides, semi-volatile organic compounds and other target analytes

Learn more at cem.com/edge

CEM Advancing Science

cem.com

© 2026 CEM Corporation. All Rights Reserved.